

For Research Use Only

EVI5L Polyclonal antibody, PBS Only

Catalog Number: 32934-1-PBS



Basic Information

Catalog Number:

32934-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG35996

GenBank Accession Number:

BC014111

GeneID (NCBI):

115704

UNIPROT ID:

Q96CN4

Full Name:

ecotropic viral integration site 5-like

Observed MW:

91,92 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

Evi5l (Ecotropic Viral Integration Site 5-Like), encodes a RabGAP that is known to block the formation of cilia when overexpressed (PMID: 26136363). EVI5L has ubiquitously low expression but is slightly higher in the testis and fetal brain.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32934-1-AP (EVI5L antibody) at dilution of 1:6000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 32934-1-PBS in a different storage buffer formulation.