

For Research Use Only

ALS2CL Polyclonal antibody, PBS Only

Catalog Number: 32937-1-PBS



Basic Information

Catalog Number:

32937-1-PBS

Size:

100ug, Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36619

GenBank Accession Number:

BC075825

GeneID (NCBI):

259173

UNIPROT ID:

Q60I27

Full Name:

ALS2 C-terminal like

Observed MW:

27kDa, 90-100kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32937-1-AP (ALS2CL antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 32937-1-PBS in a different storage buffer formulation.