

For Research Use Only

RNF185 Polyclonal antibody, PBS Only

Catalog Number: 32956-1-PBS



Basic Information

Catalog Number:

32956-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG38937

GenBank Accession Number:

NM_152267.3

GeneID (NCBI):

91445

UNIPROT ID:

Q96GF1

Full Name:

ring finger protein 185

Calculated MW:

20kDa 192aa

Observed MW:

22 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

RNF185 is a RING finger domain-containing ubiquitin ligase implicated in ER-associated degradation. It is the MOM E3 ligase responsible for regulation of mitochondrial autophagy (PMID: 21931693).

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32956-1-AP (RNF185 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 32956-1-PBS in a different storage buffer formulation.