

For Research Use Only

ADAMTSL5 Polyclonal antibody

Catalog Number:32965-1-AP



Basic Information

Catalog Number:

32965-1-AP

Size:

150ul , Concentration: 300 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG37537

GenBank Accession Number:

NM_213604

GeneID (NCBI):

339366

UNIPROT ID:

Q6ZMM2

Full Name:

ADAMTS-like 5

Calculated MW:

53 kDa

Observed MW:

75 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : COLO 320 cells, Caco-2 cells, HeLa cells, mouse skeletal muscle tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32965-1-AP (ADAMTSL5 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.