

For Research Use Only

B3GALT5 Polyclonal antibody, PBS Only

Catalog Number: 33022-1-PBS



Basic Information

Catalog Number: 33022-1-PBS	GenBank Accession Number: BC104862	Purification Method: Antigen affinity Purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 10317	
Source: Rabbit	UNIPROT ID: Q9Y2C3	
Isotype: IgG	Full Name: UDP-Gal:betaGlcNAc beta 1,3-galactosyltransferase, polypeptide 5	
Immunogen Catalog Number: AG37977	Calculated MW: 310 aa, 36 kDa	
	Observed MW: 44 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human, mouse, rat

Background Information

B3GALT5 (β -1,3-galactosyltransferase 5) is a key glycosyltransferase enzyme. It is responsible for catalyzing the synthesis of the terminal glycan structure, the sialyl Lewis X antigen. This antigen serves as an important ligand for selectin family proteins, mediating cell adhesion in both physiological and pathological processes. Its core function lies in regulating cell migration. In immunology, it influences lymphocyte homing; in oncology, its aberrant upregulation promotes the synthesis of sLe^x, leading to the adhesion of cancer cells to vascular endothelium via the "selectin-ligand" recognition system. This represents a critical initial step in tumor hematogenous metastasis, such as to the liver or lungs.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33022-1-AP (B3GALT5 antibody) at dilution of 1:15000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 33022-1-PBS in a different storage buffer formulation.