

For Research Use Only

FAM166A Polyclonal antibody

Catalog Number: 33044-1-AP



Basic Information

Catalog Number:

33044-1-AP

Size:

150ul , Concentration: 400 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG37804

GenBank Accession Number:

NM_001001710

GeneID (NCBI):

401565

UNIPROT ID:

Q6J272

Full Name:

family with sequence similarity 166, member A

Calculated MW:

317aa, 36 kDa

Observed MW:

36 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:1000-1:6000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse testis tissue,

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Mouse testis lysate was subjected to SDS PAGE followed by western blot with 33044-1-AP (FAM166A antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.