

For Research Use Only

FAM166A Polyclonal antibody, PBS Only

Catalog Number:33044-1-PBS



Basic Information

Catalog Number:

33044-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG37804

GenBank Accession Number:

NM_001001710

GeneID (NCBI):

401565

UNIPROT ID:

Q6J272

Full Name:

family with sequence similarity 166, member A

Calculated MW:

317aa, 36 kDa

Observed MW:

36 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Mouse testis lysate was subjected to SDS PAGE followed by western blot with 33044-1-AP (FAM166A antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 33044-1-PBS in a different storage buffer formulation.