

For Research Use Only

BXDC5 Polyclonal antibody

Catalog Number: 33133-1-AP



Basic Information

Catalog Number: 33133-1-AP	GenBank Accession Number: BC016051	Purification Method: Antigen affinity Purification
Size: 150ul , Concentration: 300 ug/ml by Nanodrop;	GeneID (NCBI): 80135	Recommended Dilutions: WB: 1:5000-1:50000 IF/ICC: 1:400-1:1600
Source: Rabbit	UNIPROT ID: Q9H9Y2	
Isotype: IgG	Full Name: brix domain containing 5	
Immunogen Catalog Number: AG37505	Calculated MW: 349 aa, 40 kDa	
	Observed MW: 48 kDa	

Applications

Tested Applications: WB, IF/ICC, ELISA	Positive Controls: WB : HeLa cells, U-251 cells, RKO cells, RT-4 cells IF/ICC : U2OS cells,
Species Specificity: human	

Background Information

BXDC5 (Brix domain-containing protein 5), also known as ribosome biogenesis protein RPF1, is a protein containing a Brix domain and is localized in the nucleolus. The Brix domain is a region that shares homology with the yeast protein Pitx1 (ribosome biogenesis protein BRX1). BXDC5 is a key factor in ribosome biogenesis, participating in the maturation of the 60S ribosomal subunit. This protein exerts its function by binding to specific rRNA structures and acts as a downstream effector molecule of the Ran signaling pathway, involved in ribosome biogenesis. In colorectal cancer (CRC), BXDC5 influences the occurrence and development of tumors by regulating ribosome biogenesis.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33133-1-AP (BXDC5 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (4% PFA) fixed U2OS cells using BXDC5 antibody (33133-1-AP) at dilution of 1:800 and CoraLite@488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red).