

For Research Use Only

# GTF3A Polyclonal antibody

Catalog Number: 33192-1-AP



## Basic Information

|   |  |  |
|---|--|--|
| <b>Catalog Number:</b><br>33192-1-AP                          | <b>GenBank Accession Number:</b><br>NM_002097          | <b>Purification Method:</b><br>Antigen affinity Purification |
| <b>Size:</b><br>150ul , Concentration: 450 ug/ml by Nanodrop; | <b>GeneID (NCBI):</b><br>2971                          | <b>Recommended Dilutions:</b><br>WB: 1:1000-1:4000           |
| <b>Source:</b><br>Rabbit                                      | <b>UNIPROT ID:</b><br>Q92664                           |  |
| <b>Isotype:</b><br>IgG  | <b>Full Name:</b><br>general transcription factor IIIA |  |
| <b>Immunogen Catalog Number:</b><br>AG35931                   | <b>Calculated MW:</b><br>41kDa                         |  |
|   | <b>Observed MW:</b><br>42 kDa                          |  |

## Applications

|  |   |
|--|---|
| <b>Tested Applications:</b><br>WB, ELISA | <b>Positive Controls:</b><br>WB : LNCaP cells, PC-3 cells |
| <b>Species Specificity:</b><br>human     |   |

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol, pH7.3  
Aliquoting is unnecessary for -20°C storage

**\*\*\* 20ul sizes contain 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33192-1-AP (GTF3A antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.