

For Research Use Only

KCTD19 Polyclonal antibody

Catalog Number:33194-1-AP



Basic Information

Catalog Number:

33194-1-AP

Size:

150ul , Concentration: 400 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36660

GenBank Accession Number:

NM_001100915.3

GeneID (NCBI):

146212

UNIPROT ID:

Q17RG1

Full Name:

potassium channel tetramerisation domain containing 19

Calculated MW:

105 kDa 926aa

Observed MW:

105 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse testis tissue, rat testis tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33194-1-AP (KCTD19 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.