

For Research Use Only

KCTD19 Polyclonal antibody, PBS Only

Catalog Number: 33194-1-PBS



Basic Information

Catalog Number:

33194-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36660

GenBank Accession Number:

NM_001100915.3

GeneID (NCBI):

146212

UNIPROT ID:

Q17RG1

Full Name:

potassium channel tetramerisation domain containing 19

Calculated MW:

105 kDa 926aa

Observed MW:

105 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33194-1-AP (KCTD19 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 33194-1-PBS in a different storage buffer formulation.