

For Research Use Only

PIM3 Polyclonal antibody, PBS Only

Catalog Number: 33202-1-PBS



Basic Information

Catalog Number:

33202-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG39664

GenBank Accession Number:

BC141855

GeneID (NCBI):

415116

UNIPROT ID:

Q86V86

Full Name:

pim-3 oncogene

Calculated MW:

326 aa, 36 kDa

Observed MW:

40 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33202-1-AP (PIM3 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 33202-1-PBS in a different storage buffer formulation.