

For Research Use Only

FAM21A Polyclonal antibody

Catalog Number: 33208-1-AP



Basic Information

Catalog Number: 33208-1-AP	GenBank Accession Number: BC082258	Purification Method: Antigen affinity Purification
Size: 150ul , Concentration: 400 ug/ml by Nanodrop;	GeneID (NCBI): 387680	Recommended Dilutions: WB: 1:2000-1:16000
Source: Rabbit	UNIPROT ID: Q641Q2	
Isotype: IgG	Full Name: family with sequence similarity 21, member A	
Immunogen Catalog Number: AG36997	Observed MW: 170 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : Caco-2 cells,
Species Specificity: human	

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33208-1-AP (FAM21A antibody) at dilution of 1:8000 incubated at room temperature for 1.5 hours.