

For Research Use Only

FAM21A Polyclonal antibody, PBS Only

Catalog Number: 33208-1-PBS



Basic Information

Catalog Number:

33208-1-PBS

Size:

100ug, Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36997

GenBank Accession Number:

BC082258

GeneID (NCBI):

387680

UNIPROT ID:

Q641Q2

Full Name:

family with sequence similarity 21,
member A

Observed MW:

170 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33208-1-AP (FAM21A antibody) at dilution of 1:8000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 33208-1-PBS in a different storage buffer formulation.