

For Research Use Only

# ZNF551 Polyclonal antibody, PBS Only

Catalog Number: 33216-1-PBS



## Basic Information

Catalog Number:

33216-1-PBS

Size:

100ug, Concentration: 1 mg/ml by  
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG35151

GenBank Accession Number:

BC005868

GeneID (NCBI):

90233

UNIPROT ID:

Q7Z340

Full Name:

zinc finger protein 551

Calculated MW:

670 aa, 78 kDa

Observed MW:

41 kDa

Purification Method:

Antigen affinity Purification

## Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human, mouse

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33216-1-AP (ZNF551 antibody) at dilution of 1:6000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 33216-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded mouse liver tissue slide using 33216-1-AP (ZNF551 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 33216-1-PBS in a different storage buffer formulation.