

For Research Use Only

TRIM77 Polyclonal antibody

Catalog Number:33260-1-AP



Basic Information

Catalog Number:

33260-1-AP

Size:

150ul , Concentration: 350 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG37467

GenBank Accession Number:

NM_001146162

GeneID (NCBI):

390231

UNIPROT ID:

I1YAP6

Full Name:

tripartite motif-containing 77

Calculated MW:

450 aa, 52 kDa

Observed MW:

70 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse pancreas tissue, mouse testis tissue, rat pancreas tissue, rat testis tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33260-1-AP (TRIM77 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.