

For Research Use Only

OPN1LW Polyclonal antibody, PBS Only

Catalog Number: 33320-1-PBS



Basic Information

Catalog Number:

33320-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG37894

GenBank Accession Number:

BC156643

GeneID (NCBI):

5956

UNIPROT ID:

P04000

Full Name:

opsin 1 (cone pigments), long-wave-sensitive

Observed MW:

40 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

IP, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



IP result of anti-OPN1LW (IP:33320-1-AP, 4ug; Detection:33320-1-AP 1:500) with mouse retina tissue lysate 720 ug. This data was developed using the same antibody clone with 33320-1-PBS in a different storage buffer formulation.

IP result of anti-OPN1LW (IP:33320-1-AP, 4ug; Detection:33320-1-AP 1:500) with rat retina tissue lysate 400 ug. This data was developed using the same antibody clone with 33320-1-PBS in a different storage buffer formulation.