

For Research Use Only

# TMCO2 Polyclonal antibody, PBS Only

Catalog Number:33338-1-PBS



## Basic Information

Catalog Number:

33338-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG37348

GenBank Accession Number:

NM\_001008740.3

GeneID (NCBI):

127391

UNIPROT ID:

Q7Z6W1

Full Name:

transmembrane and coiled-coil domains 2

Calculated MW:

20kDa 182aa

Observed MW:

20-25 kDa

Purification Method:

Antigen affinity Purification

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33338-1-AP (TMCO2 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 33338-1-PBS in a different storage buffer formulation.