

For Research Use Only

POU6F2 Polyclonal antibody

Catalog Number:33362-1-AP



Basic Information

Catalog Number:

33362-1-AP

Size:

150ul , Concentration: 400 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG38018

GenBank Accession Number:

NM_001166018

GeneID (NCBI):

11281

UNIPROT ID:

P78424

Full Name:

POU class 6 homeobox 2

Calculated MW:

691aa, 73 kDa

Observed MW:

73 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse brain tissue, mouse cerebellum tissue, rat brain tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33362-1-AP (POU6F2 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.