

For Research Use Only

POU6F2 Polyclonal antibody, PBS Only

Catalog Number: 33362-1-PBS



Basic Information

Catalog Number:

33362-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG38018

GenBank Accession Number:

NM_001166018

GeneID (NCBI):

11281

UNIPROT ID:

P78424

Full Name:

POU class 6 homeobox 2

Calculated MW:

691aa, 73 kDa

Observed MW:

73 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33362-1-AP (POU6F2 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 33362-1-PBS in a different storage buffer formulation.