

For Research Use Only

FCRL2 Polyclonal antibody

Catalog Number: 33416-1-AP



Basic Information

Catalog Number:

33416-1-AP

Size:

150ul, Concentration: 470 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG4773

GenBank Accession Number:

NM_030764.3

GeneID (NCBI):

79368

UNIPROT ID:

Q96LA5-1

Full Name:

Fc receptor-like 2

Calculated MW:

56 kDa

Observed MW:

60-70 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HeLa cells,

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 33416-1-AP (FCRL2 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.