

For Research Use Only

VSTM2B Polyclonal antibody, PBS Only

Catalog Number:33423-1-PBS



Basic Information

Catalog Number:

33423-1-PBS

Size:

100ug , Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG6734

GenBank Accession Number:

NM_001146339.1

GeneID (NCBI):

342865

UNIPROT ID:

A6NLU5

Full Name:

V-set and transmembrane domain
containing 2B

Calculated MW:

30 kDa

Observed MW:

30 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Fetal human brain tissue were subjected to SDS PAGE followed by western blot with 33423-1-AP (VSTM2B antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 33423-1-PBS in a different storage buffer formulation.