

For Research Use Only

CCDC58 Polyclonal antibody, PBS Only

Catalog Number: 33455-1-PBS



Basic Information

Catalog Number:

33455-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG38762

GenBank Accession Number:

BC062725

GeneID (NCBI):

131076

UNIPROT ID:

Q4VC31

Full Name:

coiled-coil domain containing 58

Observed MW:

16 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, IP, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33455-1-AP (CCDC58 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 33455-1-PBS in a different storage buffer formulation.



IP result of anti-CCDC58 (IP:33455-1-AP, 4ug; Detection:33455-1-AP 1:1000) with rat kidney tissue lysate 3160 ug. This data was developed using the same antibody clone with 33455-1-PBS in a different storage buffer formulation.