

For Research Use Only

MED7 Polyclonal antibody, PBS Only

Catalog Number: 33458-1-PBS



Basic Information

Catalog Number:

33458-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG38214

GenBank Accession Number:

BC006968

GeneID (NCBI):

9443

UNIPROT ID:

O43513

Full Name:

mediator complex subunit 7

Calculated MW:

233 aa, 27 kDa

Observed MW:

30 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

MED7 is a crucial subunit of the Mediator complex, a protein complex that regulates gene transcription by interacting with RNA polymerase II.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33458-1-AP (MED7 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 33458-1-PBS in a different storage buffer formulation.