

For Research Use Only

FAM128B Polyclonal antibody

Catalog Number:33491-1-AP



Basic Information

Catalog Number: 33491-1-AP	GenBank Accession Number: BC066296	Purification Method: Antigen affinity Purification
Size: 150ul , Concentration: 450 ug/ml by Nanodrop;	GeneID (NCBI): 80097	Recommended Dilutions: WB: 1:1000-1:4000
Source: Rabbit	UNIPROT ID: Q6NZ67	
Isotype: IgG	Full Name: family with sequence similarity 128, member B	
Immunogen Catalog Number: AG38628	Observed MW: 16-20 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HeLa cells, HepG2 cells, L02 cells
Species Specificity: human	

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1%BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33491-1-AP (FAM128B antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.