

For Research Use Only

FAM128B Polyclonal antibody, PBS Only

Catalog Number: 33491-1-PBS



Basic Information

Catalog Number: 33491-1-PBS	GenBank Accession Number: BC066296	Purification Method: Antigen affinity Purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 80097	
Source: Rabbit	UNIPROT ID: Q6NZ67	
Isotype: IgG	Full Name: family with sequence similarity 128, member B	
Immunogen Catalog Number: AG38628	Observed MW: 16-20 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33491-1-AP (FAM128B antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 33491-1-PBS in a different storage buffer formulation.