

For Research Use Only

Angiotensin 4 Polyclonal antibody, PBS Only

Catalog Number:33510-1-PBS



Basic Information

Catalog Number:

33510-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36152

GenBank Accession Number:

BC111976

GeneID (NCBI):

51378

UNIPROT ID:

Q9Y264

Full Name:

angiotensin 4

Observed MW:

57 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33510-1-AP (ANGPT4 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 33510-1-PBS in a different storage buffer formulation.