

For Research Use Only

TMEM256 Polyclonal antibody, PBS Only

Catalog Number: 33511-1-PBS



Basic Information

Catalog Number: 33511-1-PBS	GenBank Accession Number: BC030270	Purification Method: Antigen affinity Purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 254863	
Source: Rabbit	UNIPROT ID: Q8N2U0	
Isotype: IgG	Full Name: chromosome 17 open reading frame 61	
Immunogen Catalog Number: AG36789	Observed MW: 12 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

TMEM256, also known as C17orf61, is a member of TMEM families, encodes a protein of ~11.7 kDa, TMEM256 is expressed in proximal renal tubular cells, and its subcellular localization is in exosomes (PMID: 39991002).

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33511-1-AP (TMEM256 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 33511-1-PBS in a different storage buffer formulation.