

For Research Use Only

SELT Polyclonal antibody

Catalog Number: 33518-1-AP



Basic Information

Catalog Number:

33518-1-AP

Size:

150ul , Concentration: 350 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG38727

GenBank Accession Number:

BC008411

GeneID (NCBI):

51714

UNIPROT ID:

P62341

Full Name:

selenoprotein T

Calculated MW:

22 kDa

Observed MW:

18 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : BxPC-3 cells, MCF-7 cells, Ramos cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33518-1-AP (SELT antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.