

For Research Use Only

S100PBP Polyclonal antibody, PBS Only

Catalog Number: 33562-1-PBS



Basic Information

Catalog Number:

33562-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG38238

GenBank Accession Number:

BC015175

GeneID (NCBI):

64766

UNIPROT ID:

Q96BU1

Full Name:

S100P binding protein

Observed MW:

55 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33562-1-AP (S100PPB antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 33562-1-PBS in a different storage buffer formulation.