

For Research Use Only

TMEM155 Polyclonal antibody, PBS Only

Catalog Number:33601-1-PBS



Basic Information

Catalog Number:

33601-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG39858

GenBank Accession Number:

NM_001384332.1

GeneID (NCBI):

132332

UNIPROT ID:

Q4W5P6

Full Name:

transmembrane protein 155

Calculated MW:

7 kDa 63aa

Observed MW:

10-15 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, IP, Indirect ELISA

Species Specificity:

human, mouse, rat, pig

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33601-1-AP (TMEM155 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 33601-1-PBS in a different storage buffer formulation.



IP result of anti-TMEM155 (IP:33601-1-AP, 4 ug; Detection:33601-1-AP 1:500) with mouse brain tissue lysate 1720 ug. This data was developed using the same antibody clone with 33601-1-PBS in a different storage buffer formulation.