

For Research Use Only

CNTN5 Polyclonal antibody, PBS Only

Catalog Number:33610-1-PBS



Basic Information

Catalog Number:

33610-1-PBS

Size:

100ug , Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG38223

GenBank Accession Number:

BC130619

GeneID (NCBI):

53942

UNIPROT ID:

O94779

Full Name:

contactin 5

Calculated MW:

121 kDa

Observed MW:

121-150 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33610-1-AP (CNTN5 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 33610-1-PBS in a different storage buffer formulation.