

For Research Use Only

TMEM260 Polyclonal antibody

Catalog Number:33682-1-AP



Basic Information

Catalog Number:

33682-1-AP

Size:

150ul , Concentration: 340 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG38950

GenBank Accession Number:

NM_017799.3

GeneID (NCBI):

54916

UNIPROT ID:

Q9NX78

Full Name:

chromosome 14 open reading frame 101

Calculated MW:

80 kDa 707aa

Observed MW:

80 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HepG2 cells, HuH-7 cells, L02 cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33682-1-AP (TMEM260 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.