

For Research Use Only

# CNTN3 Polyclonal antibody

Catalog Number:33720-1-AP



## Basic Information

**Catalog Number:**

33720-1-AP

**Size:**

150ul , Concentration: 500 ug/ml by Nanodrop;

**Source:**

Rabbit

**Isotype:**

IgG

**Immunogen Catalog Number:**

EG5928

**GenBank Accession Number:**

NM\_020872.2

**GeneID (NCBI):**

5067

**UNIPROT ID:**

Q9P232

**Full Name:**

contactin 3 (plasmacytoma associated)

**Calculated MW:**

113 kDa

**Observed MW:**

140 kDa

**Purification Method:**

Antigen affinity Purification

**Recommended Dilutions:**

WB: 1:500-1:3000

## Applications

**Tested Applications:**

WB, ELISA

**Species Specificity:**

human, mouse, rat

**Positive Controls:**

WB : mouse brain tissue, mouse cerebellum tissue, rat brain tissue, rat cerebellum tissue

## Storage

**Storage:**

Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33720-1-AP (CNTN3 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.