

For Research Use Only

ZRSR2 Polyclonal antibody

Catalog Number: 33753-1-AP



Basic Information

Catalog Number:

33753-1-AP

Size:

150ul , Concentration: 420 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36984

GenBank Accession Number:

BC050451

GeneID (NCBI):

8233

UNIPROT ID:

Q15696

Full Name:

zinc finger (CCCH type), RNA-binding motif and serine/arginine rich 2

Calculated MW:

58 kDa

Observed MW:

60-70 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HEK-293 cells, K-562 cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33753-1-AP (ZRSR2 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.