

For Research Use Only

AGK Polyclonal antibody

Catalog Number: 33759-1-AP



Basic Information

Catalog Number: 33759-1-AP	GenBank Accession Number: BC022777	Purification Method: Antigen affinity Purification
Size: 150ul , Concentration: 370 ug/ml by Nanodrop;	GeneID (NCBI): 55750	Recommended Dilutions: WB: 1:2000-1:12000 IHC: 1:50-1:500
Source: Rabbit	UNIPROT ID: Q53H12	
Isotype: IgG	Full Name: acylglycerol kinase	
Immunogen Catalog Number: AG37404	Calculated MW: 47 kDa	
	Observed MW: 47 kDa	

Applications

Tested Applications:
WB, IHC, ELISA

Species Specificity:
human, mouse, rat

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

WB : HEK-293T cells, HeLa cells, K-562 cells, L02 cells, mouse brain tissue, rat brain tissue

IHC : mouse brain tissue,

Background Information

AGK (Acylglycerol Kinase), also designated as MULK (multi-substrate lipid kinase), is a member of the AGK protein family and functions as a lipid kinase localized to the inner mitochondrial membrane. This protein exhibits a distinctive dual functionality. First, acting as a lipid kinase, it catalyzes the phosphorylation of monoacylglycerol and diacylglycerol to yield lysophosphatidic acid (LPA) and phosphatidic acid (PA), respectively, thereby participating in phospholipid metabolism and lipid-mediated signal transduction. Second, independent of its kinase activity, AGK serves as a core subunit of the TIM22 complex, where it is responsible for the translocation and insertion of multiple mitochondrial carrier proteins (e.g., ANT1) into the inner mitochondrial membrane.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33759-1-AP (AGK antibody) at dilution of 1:6000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using 33759-1-AP (AGK antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).