

For Research Use Only

C10orf35 Polyclonal antibody, PBS Only

Catalog Number:33819-1-PBS



Basic Information

Catalog Number:

33819-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36807

GenBank Accession Number:

BC013587

GeneID (NCBI):

219738

UNIPROT ID:

Q96D05

Full Name:

chromosome 10 open reading frame

35

Observed MW:

15 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33819-1-AP (C10orf35 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 33819-1-PBS in a different storage buffer formulation.