

For Research Use Only

ZNF157 Polyclonal antibody, PBS Only

Catalog Number: 33838-1-PBS



Basic Information

Catalog Number: 33838-1-PBS	GenBank Accession Number: BC075003	Purification Method: Antigen affinity Purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 7712	
Source: Rabbit	UNIPROT ID: P51786	
Isotype: IgG	Full Name: zinc finger protein 157	
Immunogen Catalog Number: AG39863	Observed MW: 55-65 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

ZNF157, also known as HZF22, is a zinc finger protein gene located at the X chromosome short arm region 11.3. As a member of the zinc finger protein family, ZNF157 primarily functions through its C2H2-type zinc finger domain. This domain enables ZNF157 to bind to DNA, thereby playing a key role in the regulation of gene expression.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33838-1-AP (ZNF157 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 33838-1-PBS in a different storage buffer formulation.