

For Research Use Only

SSBP3 Polyclonal antibody, PBS Only

Catalog Number: 33854-1-PBS



Basic Information

Catalog Number: 33854-1-PBS	GenBank Accession Number: BC066365	Purification Method: Antigen affinity Purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 23648	
Source: Rabbit	UNIPROT ID: Q9BWW4	
Isotype: IgG	Full Name: single stranded DNA binding protein 3	
Immunogen Catalog Number: AG40791	Calculated MW: 40 kDa	
	Observed MW: 40 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human, mouse, rat, pig

Background Information

SSBP 3 (single-stranded DNA binding protein 3) is a kind of nuclear transcription co-activator. Its core function is to regulate embryonic development, islet β cell function and gene expression by binding single-stranded DNA and participating in transcription regulation complex. Its abnormality is related to diseases such as metabolism and tumor.

Storage

Storage:
Store at -80°C .
Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33854-1-AP (SSBP3 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 33854-1-PBS in a different storage buffer formulation.