

For Research Use Only

ACBD7 Polyclonal antibody

Catalog Number:33855-1-AP



Basic Information

Catalog Number: 33855-1-AP	GenBank Accession Number: BC029526	Purification Method: Antigen affinity Purification
Size: 150ul , Concentration: 580 ug/ml by Nanodrop;	GeneID (NCBI): 414149	
Source: Rabbit	UNIPROT ID: Q8N6N7	
Isotype: IgG	Full Name: acyl-Coenzyme A binding domain containing 7	
Immunogen Catalog Number: AG40829	Calculated MW: 88 aa, 10 kDa	
	Observed MW: 10 kDa	

Applications

Tested Applications:
ELISA
Species Specificity:
human

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1%BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data