

For Research Use Only

# ACBD7 Polyclonal antibody, PBS Only

Catalog Number:33855-1-PBS



## Basic Information

Catalog Number:

33855-1-PBS

Size:

100ug , Concentration: 1 mg/ml by  
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG40829

GenBank Accession Number:

BC029526

GeneID (NCBI):

414149

UNIPROT ID:

Q8N6N7

Full Name:

acyl-Coenzyme A binding domain  
containing 7

Calculated MW:

88 aa, 10 kDa

Observed MW:

10 kDa

Purification Method:

Antigen affinity Purification

## Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data