

For Research Use Only

SFRS4 Polyclonal antibody

Catalog Number: 33904-1-AP



Basic Information

Catalog Number:

33904-1-AP

Size:

150ul , Concentration: 360 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG39084

GenBank Accession Number:

BC002781

GeneID (NCBI):

6429

UNIPROT ID:

Q08170

Full Name:

splicing factor, arginine/serine-rich 4

Calculated MW:

57 kDa

Observed MW:

50-60 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HEK-293 cells, Jurkat cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33904-1-AP (SFRS4 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.