

For Research Use Only

CENPS-CORT Polyclonal antibody

Catalog Number:33963-1-AP



Basic Information

Catalog Number: 33963-1-AP	GenBank Accession Number: GeneID (NCBI): 100526739	Purification Method: Antigen affinity Purification
Size: 150ul , Concentration: 360 ug/ml by Nanodrop;	UNIPROT ID: Q8N2Z9	Recommended Dilutions: WB: 1:1000-1:4000
Source: Rabbit	Full Name: CENPS-CORT readthrough	
Isotype: IgG	Calculated MW: 15kDa	
Immunogen Catalog Number: AG40952	Observed MW: 19 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : Jurkat cells, K-562 cells, THP-1 cells
Species Specificity: human	

Background Information

CENPS-CORT is a centromeric protein complex composed of core subunits such as CENPS and CENPT, playing a critical role in chromosome biology. This complex localizes to the centromere region and directly participates in the assembly and maintenance of the centromere-kinetochore structure, ensuring accurate chromosome segregation during both mitosis and meiosis. Its functions primarily include promoting nucleosome assembly at the centromere, maintaining the three-dimensional structural integrity of the centromere, and ensuring the fidelity of chromosome separation by regulating microtubule-kinetochore attachments. Studies have shown that dysfunction of the CENPS-CORT complex disrupts the accuracy of chromosome segregation, leading to chromosomal numerical abnormalities such as aneuploidy, which is closely associated with tumorigenesis and reproductive disorders. It is therefore an important subject in research on cell cycle regulation and genomic stability.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33963-1-AP (CENPS-CORT antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.