

For Research Use Only

ST6GAL1 Polyclonal antibody

Catalog Number: 34030-1-AP



Basic Information

Catalog Number: 34030-1-AP	GenBank Accession Number: NM_001113344.1	Purification Method: Antigen affinity Purification
Size: 150ul , Concentration: 420 ug/ml by Nanodrop;	GeneID (NCBI): 25197	Recommended Dilutions: WB: 1:500-1:2000
Source: Rabbit	UNIPROT ID: P13721-1	
Isotype: IgG	Full Name: ST6 beta-galactosamide alpha-2,6-sialyltransferase 1	
Immunogen Catalog Number: EG5252	Calculated MW: 47 kDa	
	Observed MW: 43-45 kDa, 50-70 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HCT 116 cells, mouse kidney tissue, rat kidney tissue
Species Specificity: human, mouse, rat	

Background Information

ST6GAL1 (β -galactoside α -2-6 sialyl transferase1; also known as ST6N or CD75) is a sialyltransferase mediating the glycosylation of proteins and lipids to form functionally important glycoproteins and glycolipids in the Golgi compartment. It is principally expressed in liver, placenta, and skeletal muscle. ST6GAL1 undergoes proteolytic process to generate soluble form from membrane form. Higher molecular weight of bands around 50-70 kDa can also be observed with glycosylation modification. (PMID: 15049997, 23358684)

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 34030-1-AP (ST6GAL1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.