

For Research Use Only

FLYWCH1 Polyclonal antibody, PBS Only

Catalog Number: 34047-1-PBS



Basic Information

Catalog Number:

34047-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG41111

GenBank Accession Number:

BC018712

GeneID (NCBI):

84256

UNIPROT ID:

Q4VC44

Full Name:

FLYWCH-type zinc finger 1

Observed MW:

85 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 34047-1-AP (FLYWCH1 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 34047-1-PBS in a different storage buffer formulation.