

For Research Use Only

KLF11 Polyclonal antibody, PBS Only

Catalog Number: 34052-1-PBS



Basic Information

Catalog Number:

34052-1-PBS

Size:

100ug, Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG40584

GenBank Accession Number:

BC063286

GeneID (NCBI):

8462

UNIPROT ID:

O14901

Full Name:

Kruppel-like factor 11

Calculated MW:

55 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human, mouse

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data