

For Research Use Only

PAXIP1 Polyclonal antibody

Catalog Number: 34133-1-AP



Basic Information

Catalog Number:

34133-1-AP

Size:

150ul , Concentration: 430 µg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG41159

GenBank Accession Number:

BC033781

GeneID (NCBI):

22976

UNIPROT ID:

Q6ZW49

Full Name:

PAX interacting (with transcription-activation domain) protein 1

Calculated MW:

118 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data