

For Research Use Only

# ROPN1B Polyclonal antibody

Catalog Number: 34136-1-AP



## Basic Information

Catalog Number:

34136-1-AP

Size:

150ul , Concentration: 450 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10787

GenBank Accession Number:

BC015413

GeneID (NCBI):

152015

UNIPROT ID:

Q9BZX4

Full Name:

ropporin, rhophilin associated protein 1B

Calculated MW:

120 aa, 24 kDa

Observed MW:

22 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:2000-1:12000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : human testis tissue, mouse testis tissue, rat testis tissue

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 34136-1-AP (RPN1B antibody) at dilution of 1:6000 incubated at room temperature for 1.5 hours.