

For Research Use Only

HEATR9 Polyclonal antibody, PBS Only

Catalog Number: 34152-1-PBS



Basic Information

Catalog Number:

34152-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG40489

GenBank Accession Number:

BC033734

GeneID (NCBI):

256957

UNIPROT ID:

A2R7Y3

Full Name:

chromosome 17 open reading frame 66

Observed MW:

70 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

HEATR9 (HEAT Repeat Containing 9, also known as C17orf66 or FLJ32830) is a protein containing HEAT repeats. It belongs to the HEAT repeat protein family, and its characteristic HEAT domain is composed of approximately 37-47 amino acid tandem repeats of amphipathic α -helical structural units, which mediate protein-protein interactions.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 34152-1-AP (C17orf66 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 34152-1-PBS in a different storage buffer formulation.