

For Research Use Only

HMX3 Polyclonal antibody

Catalog Number: 34183-1-AP



Basic Information

Catalog Number: 34183-1-AP	GenBank Accession Number: NM_001105574	Purification Method: Antigen affinity Purification
Size: 150ul , Concentration: 460 ug/ml by Nanodrop;	GeneID (NCBI): 340784	Recommended Dilutions: WB: 1:2000-1:10000 IHC: 1:200-1:800 IF/ICC: 1:300-1:1200
Source: Rabbit	UNIPROT ID: A6NHT5	
Isotype: IgG	Full Name: H6 family homeobox 3	
Immunogen Catalog Number: AG41398	Calculated MW: 357aa, 38 kDa	
	Observed MW: 38 kDa	

Applications

Tested Applications: WB, IHC, IF/ICC, ELISA	Positive Controls: WB : mouse liver tissue, rat liver tissue IHC : mouse liver tissue, IF/ICC : HeLa cells,
Species Specificity: human, mouse, rat	
Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0	

Background Information

H6 family homeobox 3 (HMX3) is also named Homeobox protein H6 family member 3 and NKX-5.1. HMX3 is identified as a protein-coding gene, and diseases associated with HMX3 primarily include upper semicircular canal dehiscence and semicircular canal dehiscence syndrome (PMID: 22779713). In addition, as a transcription factor, HMX3 is involved in the regulation of nerve cell types, which is essential for the development of the inner ear and hypothalamus (PMID: 27151948, PMID: 15363417). HMX3 is lowly expressed in CRC, and HMX3 expression was negatively correlated with poor prognosis of CRC patients (PMID: 36204976).

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 34183-1-AP (HMx3 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded mouse liver tissue slide using 34183-1-AP (HMx3 antibody) at dilution of 1:400 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using HMx3 antibody (34183-1-AP) at dilution of 1:600 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-phalloidin (red).